Products	Prime Directions	Directions for Use	Cautions
<p>Nexus TKO®-6P</p> <p>Nexus NIS®-6P</p> <p>Extension sets</p> <p><i>Note: For sets with filters please see below for specifics on filter</i></p>	<ul style="list-style-type: none"> Remove device / set from packaging using aseptic technique Swab each device in a circular motion for minimum of 10 sec with alcohol prep pad and allow to dry When a TKO or NIS is present: <ul style="list-style-type: none"> Completely prime the device / set with a flushing syringe by tilting the flushing syringe up to first expel the air from the syringe and continue to purge air through the device / set. Do not use an administration set to prime the device / set. When connecting directly to luer connector: <ul style="list-style-type: none"> Completely prime the set by connecting to primed IV administration set or using syringe. Repeat for each device / line as needed 	<ul style="list-style-type: none"> Remove protective cap, then push and twist the male luer onto the patient access device to secure connection Prior to access, always swab device as described Attach a primed IV set or syringe to the device For blood draws, attach a blood collection tube or syringe to device after disinfection and aspirate Flush the device and/or set after each use with flushing syringe Close the clamp (if present) on the patient access device then disconnect 	<ul style="list-style-type: none"> Device may not be bonded to the set, check labeling. If using a non-bonded set, tighten the device connection by firmly twisting clockwise prior to use. After blood draws, use flushing syringes as stated in directions. Do not use administration set for flushing procedure. Do not use needles, metal/blunt cannulas, non ISO compliant luers or luer connections with visible defects when accessing the device. Doing so may result in leakage and/or failure of the device. Slip luer connections should not be left unattended due to potential for disconnection. Do not use TKO device for direct pressure monitoring applications. For power injections (only if product is indicated as pressure-rated on unit label): <ul style="list-style-type: none"> Ensure catheter patency prior to power injection Verify power injector setting does not exceed the maximum pressure rating of the catheter or tubing sets
<p>Nexus TKO® One-Way Check Valve</p> <p>Extension sets</p>	<ul style="list-style-type: none"> Remove device / set from packaging using aseptic technique Swab each device in a circular motion for minimum of 10 sec with alcohol prep pad and allow to dry Completely prime the device / set with a flushing syringe by tilting the flushing syringe up to first expel the air from the syringe and continue to purge air through the device / set. Do not use an administration set to prime the device / set. Repeat for each device / line as needed 	<ul style="list-style-type: none"> Remove protective cap, then push and twist the male luer onto the patient access device to secure connection Prior to access, always swab device as described Attach a primed IV set or syringe to the device Flush the device and/or set after each use with flushing syringe Close the clamp (if present) on the patient access device then disconnect 	<ul style="list-style-type: none"> One-way check valves not intended for gravity infusion. Do not use needles, metal/blunt cannulas, non ISO compliant luers or luer connections with visible defects when accessing the device. Doing so may result in leakage and/or failure of the device. Slip luer connections should not be left unattended due to potential for disconnection.
<p>Nexus TKO® and One-Way Check Valve manifold sets</p>	<ul style="list-style-type: none"> Remove set from packaging using aseptic technique Swab each device in a circular motion for minimum of 10 sec with alcohol prep pad and allow to dry when applicable Prime each port with flushing syringe by tilting the flushing syringe up to first expel the air from the syringe and continue to purge air through the port Remove the protective cap from the female luer and prime the set with a flushing syringe or fully primed IV administration set 	<ul style="list-style-type: none"> Remove protective cap, then push and twist the male luer onto the patient access device to secure connection Prior to access, always swab device as described Attach a primed IV set to port Flush the device and/or set after each use with flushing syringe Close the clamp (if present) on the patient access device then disconnect 	<ul style="list-style-type: none"> One-way check valves not intended for gravity infusion. Do not use needles, metal/blunt cannulas, non ISO compliant luers or luer connections with visible defects when accessing the device. Doing so may result in leakage and/or failure of the device. Slip luer connections should not be left unattended due to potential for disconnection. Ensure open ports are capped when not in use. Do not use TKO device for direct pressure monitoring applications.
<p>Extension sets with filters</p>	<ul style="list-style-type: none"> Remove set from packaging using aseptic technique Swab each device in a circular motion for minimum of 10 sec with alcohol prep pad and allow to dry When a TKO or NIS is present: <ul style="list-style-type: none"> Completely prime the device / set with a flushing syringe by tilting the flushing syringe up to first expel the air from the syringe and continue to purge air through the set. Tap the filter to dislodge any air bubbles adhering to the membrane. When connecting directly to luer connector: <ul style="list-style-type: none"> Completely prime the set by connecting to primed IV administration set. Invert and tap the filter or orient such that inlet of the filter is vertically above the outlet to expel the air. Repeat for each device / line as needed 	<ul style="list-style-type: none"> Remove protective cap, then push and twist the male luer onto the patient access device to secure connection Prior to access, always swab device as described Attach a primed IV set or syringe to the device Flush the device and/or set after each use with flushing syringe Close the clamp (if present) on the patient access device then disconnect 	<ul style="list-style-type: none"> Filter will not pass whole blood, plasma products or medications not totally dissolved in solution. Do not subject filters or housings to pressures greater than 45 psi. Do not use 0.22 micron filters for lipids. Attempting to clear a filter with a syringe may lead to vent leakage. See above for additional cautions on sets containing the TKO/NIS devices.


General Cautions:

- Do not allow air to be trapped in set where applicable
- Sterile, non-pyrogenic in unopened, undamaged package
- Discard set if protective cap(s) is not in place
- Single use, do not re-sterilize
- Replace device / set in accordance to CDC guidelines or other applicable facility protocol
- Rx only

General Notes:

- Extension sets with non-bonded devices are noted in the description of the product.
- Extension sets intended for pressure infusion are noted as pressure-rated in the description of the product. Pressure infusion is rated to 325 psi or 10mL sec maximum flow rate.
- Natural rubber latex is not used in the manufacturing of the products.
- Products are DEHP-Free
- For patent information, see www.nexusmedical.com/patents
- NEXUS TKO®, NEXUS NIS®, IV Plus®, PRODUCTIVE FLUSH®, FIREWALL FOR YOUR CATHETER®, PIV5R®, PIV FIVE RIGHTS®, & RIGHT APPROACH. RIGHT RESULTS® are registered trademarks of Nexus Medical LLC

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